

### Measles IgG

**Effective Date: September 1, 2020**

The Special Chemistry section of Beaumont Laboratory Royal Oak currently performs testing for the qualitative detection of IgG antibody to the measles virus using an automated immunoassay method on the DiaSorin Liaison platform. Commencing September 1, 2020, new cut-off values that confer an increase in test sensitivity will be adopted – see below.

#### Measles IgG Cut-off values

- Units: AU (arbitrary units) / mL

Result	NEW Cut-Off Values	OLD Cut-Off Values
Negative	< 13.5	< 25.0
Equivocal	13.5 – 16.4	25.0 – 29.9
Positive	> 16.4	> 29.9

#### Measles IgG Ordering Information:

Epic Test Code: LAB6898

Specimen Collection Requirement: One gold-top Serum Separator Tube (SST)

*If you have questions, please contact Client Services (1-800-551-0488, option 5).*

For additional information, access the Laboratory Test Directory at:

<http://beaumontlaboratory.com/test-lab-directory>

---

**Date submitted:** August 28, 2020

**Submitted By:** Gabriel Maine, PhD – Technical Director  
Elizabeth Sykes, MD – Medical Director

Department of Clinical Pathology  
Section of Special Chemistry, Royal Oak

---

*Beaumont Laboratory*  
Customer Service  
1-800-551-0488