

Green-top (lithium heparin additive) Microtainer® with Plasma Separator Now Available

Effective Date: July 1, 2021

Effective July 2021, a new blood collection tube type – green-top (lithium heparin additive) Microtainer® with plasma separator will be available for pediatric, NICU and other units at Farmington Hills, Royal Oak, Grosse Pointe, and Troy. The PeopleSoft ordering number for this new tube type is 2068341.

Please note that this will be not an acceptable tube type for Cardiac Troponin (EPIC: LAB7028, TROPG).



If you have questions, please contact Client Services (1-800-551-0488, option 5).

Laboratory Test Directory: <http://beaumontlaboratory.com/test-lab-directory>.

Date submitted: June 29, 2021

Submitted by: Katie (Qian) Sun, PhD, Technical Director, Automated Chemistry Core Lab, RO
Steven Truscott, PhD, Technical Director, Special Chemistry, STAT Lab,
Point of Care Testing
Elizabeth Sykes, MD, System Medical Director, Clinical Chemistry

Beaumont Laboratory

28050 Grand River Ave.
Farmington Hills, MI 48336

468 Cadieux Road
Grosse Pointe, MI 48230

3601 West 13 Mile Road
Royal Oak, MI 48073

44201 Dequindre Road
Troy, MI 48085

www.beaumontlaboratory.com